

	HAL Lifesciences	Date: 16/12/2015	
		Doc no.:	

Spike

Myfillphase

Functional Design Specification for Manufacturing
Control System

	HAL Lifesciences	Date: 16/12/2015	
		Doc no.:	


APPROVALS

Functional Area	Required Approver	Date
Author / Originator	Print Name:	
	Sign Name:	
Engineering	Print Name:	
	Sign Name:	
Operations	Print Name:	
	Sign Name:	
Automation	Print Name:	
	Sign Name:	
QAV	Print Name:	
	Sign Name:	

	HAL Lifesciences	Date: 16/12/2015	
		Doc no.:	

Revision History

Revision	Date	Author	Description

	<p style="text-align: center;">HAL Lifesciences</p>	Date: 16/12/2015	
		Doc no.:	

This phase template is for filling a vessel. It has built in failure monitor logic and watchdog timers. It requires the following in order to successfully run:

unit with specific parameters as outlined below

creation of CM and EM physical model as outlined below

creation of formula with specific parameters.

Finally, the phase must be encased in a procedure to be executed via the campaign manager.

Table of Contents

APPROVALS	2
Revision History	3
Table of Contents.....	4
Reference Documentation.....	7
1. Myfillphase.....	8
1.1. Held	8
1.1.1. S0 Logic.....	8
1.1.2. T0 Logic	8
1.1.3. T1 Logic	9
1.2. Restarting	9
1.2.1. S0 Logic.....	9
1.2.2. T0 Logic	10
1.2.3. T1 Logic	10
1.3. Running	11



HAL Lifesciences

Date: 16/12/2015

Doc no.:

1.3.1.	S0 Logic.....	11
1.3.2.	S1 Logic.....	13
1.3.3.	S2 Logic.....	13
1.3.4.	S3 Logic.....	14
1.3.5.	S4 Logic.....	14
1.3.6.	S5 Logic.....	15
1.3.7.	S6 Logic.....	15
1.3.8.	T0 Logic	16
1.3.9.	T1 Logic	16
1.3.10.	T2 Logic	16
1.3.11.	T3 Logic	16
1.3.12.	T4 Logic	17
1.3.13.	T5 Logic	17
1.3.14.	T6 Logic	17
1.3.15.	T7 Logic	17
1.3.16.	Failure Monitor Logic.....	18
2.	Formula Name.....	19

	HAL Lifesciences	Date: 16/12/2015	
		Doc no.:	


	HAL Lifesciences	Date: 16/12/2015	
		Doc no.:	

Reference Documentation

Document	Ref.	Title

1. Myfillphase

1.1.1. Held



1.1.1.1. S0 Logic

```

#####
# Step: Held/S0
#####

def main():
    """Main Function"""

    # failure monitor timer hold.

    # get a handle to the timer
    mytimer = unit.SupportModule.GetTimer('FailureMonitorTimer')

    # if in HELD state, stop timer
    mytimer.Stop()

    # print status to consle window

    info('phase HELD')
  
```

1.1.1.2. T0 Logic

```

#####
  
```




```

# Transition: Held/T0
#####

def main():
    """Main Function"""

    # Transition code...

```

1.1.1.3. T1 Logic

```


#####
# Transition: Held/T1
#####

def main():
    """Main Function"""

    # Transition code...

```

1.1.2. Restarting



1.1.2.1. S0 Logic

```

#####
# Step: Restarting/S0
#####

def main():
    """Main Function"""

```

	HAL Lifesciences	Date: 16/12/2015	
		Doc no.:	

```

# get a handle to the failure monitor timer

mytimer = unit.SupportModule.GetTimer('FailureMonitorTimer')

# reset the failure monitor watchdog timer upon restarting the phase
mytimer.Reset()

# start the timer again
mytimer.Start()

```

1.1.2.2. T0 Logic

```

#####
# Transition: Restarting/T0
#####

def main():
    """Main Function"""

    # Transition code...

```

1.1.2.3. T1 Logic

```


#####
# Transition: Restarting/T1
#####

def main():
    """Main Function"""

    # Transition code...

```

1.1.3. Running



1.1.3.1. S0 Logic

#####



HAL Lifesciences

Date: 16/12/2015

Doc no.:

```
# Step: 7/S0
#####

def main():
    """Main Function"""

    # There are a number of requirements that must be met in order to run this phase:
    #
    # 1. It must run on a unit with the following unit parameters:
    # em_vent ;          this is the name of a tank vent equipment module (string)
    # em_tank_fill ;    this is the name of a tank fill equipment module (string)
    # em_water_drop ;   this is the name of a water supply equipment module that draws
from a water utility distribution header.(string)
    # leveltransmitter; this is the analogue input instrument that reads the level in
the tank. (string)
    # watchdogtimer ;   this is the maximum duration that a fill should take for this
vessel. Phase goes to HELD state, if fill setpoint not achieved in time. (double - secs)
    #
    # 2.In addition, the unit must encompass a number of equipment modules .These
modules are identified through the unit parameters or formula parameters;
    # Here, we have listed out the equipment module names directly in the unit module.
    #
    # 3.The following equipment module types and equipment module methods are required
(if you dont have them already, you will have to create them):
    #
    # em_vent (with 'isolate', 'vent', 'pressurise' methods)
    # em_tank_fill (with 'isolate', 'fill', 'cip' methods)
    # em_water_drop (with 'isolate', 'fill', 'flush' methods)
    #
    # 4. A Formula containing the following parameters must be within campaign scope ;
for example the formula must be associated with this phase or one of its parents.
    #
    # TankFillSetpoint - Fill the tank to this setpoint.(process parameter type, double
value type, units perC)
    # FillWatchdogtimer - For the failure monitor logic; if the tank is not filled
within a particular time, then raise error and drive phase to hold. (step parameter type,
double value type, units secs)
    #
    # 5. The watchdog timer value (in the formula), must be written to a unit parameter
in order for the failure monitor to read it.

    #get formula watchdogtimer and store to unitwatchdog timer

    myformulawatchdogtimer = getStepParameter('FillWatchdogtimer')

    unit.SupportModule.UpdateParameter('watchdogtimer', myformulawatchdogtimer.Value)
```



1.1.3.2. S1 Logic

```
#####  
# Step: 7/S1  
#####  
  
def main():  
    """Main Function"""  
  
    #get fill setpoint from formula  
  
    fillsetpoint = getProcessParameter('TankFillSetpoint')  
  
    # write to console window  
  
    info('about to start filling tank to setpoint: ' + str(fillsetpoint.Value))
```

1.1.3.3. S2 Logic

```
#####  
# Step: Running/S2  
#####  
  
def main():  
    """Main Function"""  
  
    # get physical model from unit parameters...  
  
    EMTankInletinstance = unit.SupportModule.PullParameter('em_tank_fill')  
    EMProcessAirinstance = unit.SupportModule.PullParameter('em_vent')  
    EMWaterDropinstance = unit.SupportModule.PullParameter('em_water_drop')  
    unitleveltransmitter = unit.SupportModule.PullParameter('leveltransmitter')  
  
    # Confirm user updated Unit parameters with the names of assigned equipment modules  
    if EMTankInletinstance.Value == 'placeholder':  
        raiseException('Unit parameters not updated from default template values.  
Please assign physical model to this unit and update unit parameters with relevant  
equipment module names.')
```



HAL Lifesciences

Date: 16/12/2015

Doc no.:

```
# sleep command is for demonstration purposes only ..
sleep(1000)

opc(EMWaterDropinstance.Value).GetMethod('flush').Run()
# sleep command is for demonstration purposes only ..
sleep(1000)

opc(EMWaterDropinstance.Value).GetMethod('fill').Run()
# sleep command is for demonstration purposes only ..
sleep(1000)

#wait for level in tank to reach setpoint

# tank level setpoint from formula
tanklevelsetpoint = getProcessParameter('TankFillSetpoint')

#use a helper function to wait until tank has reached setpoint

waitfor(untileveltransitter.Value, 'Value', Condition.GreaterThanOrEqual,
tanklevelsetpoint.Value)
```

1.1.3.4. S3 Logic

```
#####
# Step: Running/S3
#####

def main():
    """Main Function"""

    # write to console window

    info('Tank level has reached setpoint')
```

1.1.3.5. S4 Logic

```
#####
# Step: Running/S4
#####

def main():
    """Main Function"""

    # get physical model from unit parameters...
```



HAL Lifesciences

Date: 16/12/2015

Doc no.:

```
EMTankInletinstance = unit.SupportModule.PullParameter('em_tank_fill')
EMProcessAirinstance = unit.SupportModule.PullParameter('em_vent')
EMWaterDropinstance = unit.SupportModule.PullParameter('em_water_drop')
unitleveltransitter = unit.SupportModule.PullParameter('leveltransmitter')

# setup physical model for tank isolation

opc(EMWaterDropinstance.Value).GetMethod('isolate').Run()
# sleep command is for demonstration purposes only ..
sleep(1000)

opc(EMTankInletinstance.Value).GetMethod('isolate').Run()
# sleep command is for demonstration purposes only ..
sleep(1000)

opc(EMProcessAirinstance.Value).GetMethod('isolate').Run()
# sleep command is for demonstration purposes only ..
sleep(1000)
```

1.1.3.6. S5 Logic

```
#####
# Step: Running/S5
#####

def main():
    """Main Function"""

    # get handle to physical model...
    EMTankInletinstance = unit.SupportModule.PullParameter('em_tank_fill')
    EMProcessAirinstance = unit.SupportModule.PullParameter('em_vent')
    EMWaterDropinstance = unit.SupportModule.PullParameter('em_water_drop')
    unitleveltransitter = unit.SupportModule.PullParameter('leveltransmitter')

    opc(EMTankInletinstance.Value).OwnerId = ''
    opc(EMProcessAirinstance.Value).OwnerId = ''
    opc(EMWaterDropinstance.Value).OwnerId = ''
    opc(unitleveltransitter.Value).OwnerId = ''
```

1.1.3.7. S6 Logic

```
#####
# Step: Running/S6
#####

def main():
    """Main Function"""

    # get handle to physical model...
    EMTankInletinstance = unit.SupportModule.PullParameter('em_tank_fill')
    EMProcessAirinstance = unit.SupportModule.PullParameter('em_vent')
```



HAL Lifesciences

Date: 16/12/2015

Doc no.:

```
EMWaterDropinstance = unit.SupportModule.PullParameter('em_water_drop')
unitleveltransmitter = unit.SupportModule.PullParameter('leveltransmitter')

opc(EMTankInletinstance.Value).OwnerId = str(this)
opc(EMProcessAirinstance.Value).OwnerId = str(this)
opc(EMWaterDropinstance.Value).OwnerId = str(this)
opc(unitleveltransmitter.Value).OwnerId = str(this)
```

1.1.3.8. T0 Logic

```
#####
# Transition: 7/T0
#####

def main():
    """Main Function"""

    # Transition code...
```

1.1.3.9. T1 Logic

```
#####
# Transition: 7/T1
#####

def main():
    """Main Function"""

    # Transition code...
```

1.1.3.10. T2 Logic

```
#####
# Transition: Running/T2
#####

def main():
    """Main Function"""

    # Transition code...
```

1.1.3.11. T3 Logic

```
#####
# Transition: Running/T3
#####

def main():
```




```
"""Main Function"""  
  
# Transition code...
```

1.1.3.12. T4 Logic

```
#####  
# Transition: Running/T4  
#####  
  
def main():  
    """Main Function"""  
  
    # Transition code...
```

1.1.3.13. T5 Logic

```
#####  
# Transition: Running/T5  
#####  
  
def main():  
    """Main Function"""  
  
    # Transition code...
```

1.1.3.14. T6 Logic

```
#####  
# Transition: Running/T6  
#####  
  
def main():  
    """Main Function"""  
  
    # Transition code...
```

1.1.3.15. T7 Logic

```
#####  
# Transition: Running/T7  
#####  
  
def main():  
    """Main Function"""  
  
    # Transition code...
```



1.1.3.16. Failure Monitor Logic

```
#####  
# Failure Monitor (FM)  
#####  
  
def main():  
    """Main Function"""  
  
    #Failure monitor would normally monitor CM & EM alarms and watchdog timeouts.It can  
drive the phase into the HELD state based on a failure.  
    #Enter alarm monitoring logic here  
    #  
    #if OPC(module).IsInAlarm:  
    #     phase.Hold()  
    #  
    #Failure Monitor timer watchdog  
    #Fillwatchdogtimer - for the failure monitor logic ; if the tank is not filled  
within a particular time, then raise error and drive phase to hold  
    #Get watchdogtimervalue from formula  
  
    watchdogtimerduration = unit.SupportModule.PullParameter('watchdogtimer')  
  
    #Create a free running timer of 1200 seconds duration. The held logic will hold the  
timer. The restart logic will reset if phase goes to HELD  
  
    if firstScan:  
    # Create the timer once only.  
        mytimer = unit.SupportModule.NewTimer('FailureMonitorTimer', 1200)  
        mytimer.Start()  
        info('failure monitor watchdog timer started')  
  
        mytimer = unit.SupportModule.GetTimer('FailureMonitorTimer')  
  
    # Print out watchdog timer value - Helpful for troubleshooting. If this detail is  
not required and the console window is unnecessarily cluttered, then just delete.  
    if firstScan == False:  
        info('timer value:'+str(mytimer.Value))  
  
    # drive the running phase to the HELD state, if the step has been running longer  
than the watchdog timeout  
    if (mytimer.Value > watchdogtimerduration.Value) and (phase.CurrentState ==  
State.Running):  
        info('phase failure')  
        phase.Hold()
```

	HAL Lifesciences	Date: 16/12/2015	
		Doc no.:	

2. Formula Name

Formula Name
myfillmformula